

Hospira MedNet® System

Wireless 802.11a/b/g Module

Introducing the New 802.11a/b/g Wireless Module

Have you connected to the wireless solution?

Take your network connectivity to the next level with the Hospira MedNet system and the 802.11a/b/g wireless standards.

Is your data secure enough?

Hospira's compliance with IEEE and WPA/WPA-2 security standards make it possible to have an open-systems approach with simple user interfaces and the security that you need.



Wireless 802.11a/b/g module for the Plum A+® Infuser*

Specifications

Connectivity

- Wireless
 - 802.11a/b/g
 - Radio
 - Modes of Operation
 - Auto Selection (802.11a/b/g)
 - 802.11a Only
 - 802.11b Only
 - 802.11g Only
 - 802.11b/g
 - Power Save Settings
 - Full Power
 - Maximum Power Save
 - Transmit Power Settings
 - Maximum Transmit
 - Power Adjustable
 - Association Methods
 - Open
 - Shared
 - Authentication Protocols
 - EAP-TLS
 - EAP-TTLS
 - EAP-MSCHAPv2
 - EAP-MD5
 - PAP
 - CHAP
 - MSCHAP
 - MSCHAPv2
 - EAP-PEAP
 - EAP-MSCHAPv2
 - EAP-MD5
 - EAP-FAST
 - Auto PAC Provisioning
 - Manual PAC Provisioning
 - Encryption Methods
 - Static WEP w/o EAP
 - WPA-PSK w/o EAP
 - TKIP
 - AES/CCMP
 - WPA
 - TKIP
 - WPA-2
 - AES/CCMP

- Wired
 - Ethernet 10/100 10BaseT
- IP Address Modes
 - Static
 - DHCP

Configuration Tool

- HTTP/HTTPS Configuration Screens
- Administration
 - Login Name and Password Configurable
- Telnet Interface ON/OFF Configurable
- FTP Interface ON/OFF Configurable
- Hospira MedNet Server
 - Server Selectable
 - SSL Mode Selectable

Wireless LAN Device Specifications

- Standards
 - IEEE 802.11a/b/g
- Transmit Power:
 - IEEE 802.11 Standards
- Antenna
 - Integrated Antenna

Product Information

- FCC ID
 - STJ-80411396001
- IC No.
 - 5627A-80411396
- Certifications
 - FCC Part 15.247, 15.407
 - IC RSS-210, RSS-102
- Product Name
 - Hospira MedNet 802.11a/b/g Wireless Upgrade Module
- Note: For use with Hospira Plum A+ Infusion System
- Infuser Model Nos.
 - 11971-04
 - 12391-04
 - 20678-04
 - 20679-04
 - 20792-04
- Module Nos.
 - 20677-04-77
 - 20791-04-77



Advancing Wellness™

*May not be available for all models