

CLAVE[®] IV Administration System

Prepare. Transport. Deliver.



Helping Protect Patients, Caregivers
and the Environment from
Exposure to Hazardous Medications

Prepare...

This CLAVE[®] IV Administration system focuses on safety



CLAVE Universal multi-dose vial access spike with integrated 0.22 micron air vent



CLAVE multi-dose vial access spike with vent and stability skirt for use on 20 mm vials



CLAVE multi-dose vial access spike with stability skirt for use on 13 mm vials

The CLAVE[®] Drug Preparation System

- Integrated filters trap particles within the vial while equalizing pressure
- Latex-free and simple to use, with no needles or caps
- Microbiologically and mechanically closed system

Transport...

This CLAVE® IV Administration system makes transportation of hazardous medications simple, safe and secure



Spiros - A device that is closed in the unactivated state, sealing the fluid path. When attached to any female luer or valve, the fluid path is opened and allows for administration or aspiration of fluid

Spiros™ Syringe Adapter for Preparation

- Perform closed draws and transfers without the risk of drips or spills
 - Opens only when attached to a female connector, including needlefree connectors
- Safely transport syringes prepared with hazardous medication throughout a facility

Spiros Tubing Sets for Delivery

- Use on any primary, secondary or extension set in place of a standard male luer
- Prevent spills from IV sets in the event of disconnect, whether accidental or intentional
- Safely transport IV sets prepared with hazardous medications throughout a facility

Spiros Technology

- Protect the integrity of the IV fluid container, whether it is on a syringe for transfer or on the end of an IV set
 - When connected to any female luer, including needlefree connectors, the fluid path opens
 - Upon disconnect, automatically self-seals and closes the system

Deliver...

This CLAVE® IV Administration system offers flexibility, as well as safety, in the delivery of hazardous medications



Gravity Administration

- Primary and secondary gravity sets are available with the CLAVE® needlefree connector and/or Spiros™ closed male connector
 - Sets may be customized to meet individual clinician's needs
 - All sets include latex-free components and DEHP-free tubing

Security

- This CLAVE IV delivery system is specifically designed for helping safely handle and deliver hazardous IV medications
 - Sets and components work with any delivery system, regardless of the pump or tubing provider, to allow for closed, needlefree access to any solution container

For more information on Advancing Wellness™, contact your Hospira representative at 1-877-9Hospira (946-7747) or visit www.hospira.com