

Same CLAVE[®] Features, but Smaller

Same CLAVE[®] Swabable Valve Features...

- Auto shut-off upon disconnection
- Threaded lock assures secure connection
- Eliminates use of sharps with IV administration set (cannot be used with a needle)
- No caps or cannulas
- Compatible with universal luer-lock fitting

...But With Smaller Design

- Low profile to help minimize patient discomfort
- Ultra-low residual volume: 0.04 mL
- Convenient finger grip
- Flow rate of 165 mL/min: greater than an 18-G needle
- Backpressure: 75 psig
- Latex-free



The LifeShield[®]
microCLAVE[®] Connector



The LifeShield® microCLAVE® Connector

List No.	Product Description	Nominal Length	Approx. Priming Capacity (mL)	Tubing ID	Case Size
12568-01	LifeShield® microCLAVE® male adapter plug. Latex-Free. Non-DEHP.	—	—	—	100
12548-01	Microbore extension set with microCLAVE®, removable slide clamp and luer slip. Latex-Free. Non-DEHP.	7" (18 cm)	0.2	0.047"	50
12549-01	Microbore extension set with microCLAVE®, female luer lock, slide clamp, and mini-spin lock. Latex-Free. Non-DEHP.	7" (18 cm)	0.3	0.047"	50
12551-01	Macrobore extension set with microCLAVE®, female luer lock, slide clamp, and Option-Lok®. Latex-Free. Non-DEHP.	7" (18 cm)	0.7	0.090"	50
12547-01	Bifurcated microbore extension set with two removable microCLAVES®, removable slide clamps, and mini-spin lock. Latex-Free. Non-DEHP.	8" (20 cm)	0.3	0.020"	50
12546-01	Microbore extension set with luer slip Veniloop® set and microCLAVE®, Latex-Free. Non-DEHP.	4" (10 cm)	0.2	0.045"	50
12545-02	Microbore extension set with microCLAVE® spin lock T-connector, removable slide clamp, and female luer lock. Latex-Free. Non-DEHP.	5" (13 cm)	0.3	0.047"	50
12544-01	Macrobore extension set with microCLAVE® spin lock T-connector, slide clamp, and female luer lock. Latex-Free. Non-DEHP.	5" (13 cm)	0.6	0.090"	50
12550-02	Microbore extension set with microCLAVE® luer slip T-connector, removable slide clamp, and removable microCLAVE®. Latex-Free. Non-DEHP.	5" (13 cm)	0.4	0.047"	50

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Please refer to the device instructions for complete information, warnings and precautions before use.

For more information on Advancing Wellness™, contact your Hospira representative at 1-877-9Hospira (946-7747) or visit www.hospira.com