


















Symbiq™ Infusion System with Hospira MedNet® Software









Ahead of the curve: Infusion technology
designed to make patient safety easier to achieve



Symbiq™ Administration Sets

List No.	Product Description	Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
Primary I.V. Sets – Macro bore						
16001-01	LifeShield® Symbiq™ Pump Set – with convertible pin, CAIR® clamp, CLAVE® Y-site, and Option-Lok®. Latex-Free.	105" (267 cm)	15	13.0	0.090	48
						
16005-01	LifeShield® Symbiq™ Pump Set – with non-DEHP tubing, convertible pin, CAIR® clamp, CLAVE® Y-site, and Option-Lok®. Latex-Free.	105" (267 cm)	15	13.0	0.090	48
						
16008-01	LifeShield® Symbiq™ Pump Set – with convertible pin, piggyback with backcheck valve, 2 CLAVE® Y-sites, CAIR® clamp, and Option-Lok®. Latex-Free.	106" (269 cm)	15	14.0	0.090	48
						
16011-01	LifeShield® Symbiq™ Pump Nutritional Set – with convertible pin, CAIR® clamp, 1.2 micron filter, slide clamp, CLAVE® Y-site, and Option-Lok®. Latex-Free.	104" (264 cm)	15	16.0	0.090	48
						
16013-01	LifeShield® Symbiq™ Pump Set – with polyethylene-lined tubing, convertible pin, backcheck valve, 2 CLAVE® Y-sites, 0.2 micron high-pressure filter, slide clamp, and piggyback with Option-Lok®. Latex-Free.	104.5" (265 cm)	15	18.6	0.100	48
						
List No.	Product Description	Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
16019-01	LifeShield® Symbiq™ Pump Set – with non-DEHP tubing, convertible pin, piggyback with backcheck valve, 2 CLAVE® Y-sites, CAIR® clamp, and Option-Lok®. Latex-Free.	106" (269 cm)	15	14.0	0.090	48
						
16021-01	LifeShield® Symbiq™ Pump Set – with non-DEHP tubing, convertible pin, backcheck valve, 3 CLAVE® Y-sites, CAIR® clamp, slide clamp, and Option-Lok®. Latex-Free.	106" (269 cm)	15	15.0	0.090	48
						
16022-01	LifeShield® Symbiq™ Pump Set – with non-DEHP tubing, convertible pin, backcheck valve, 2 CLAVE® Y-sites, CAIR® clamp, 0.2 micron filter, slide clamp, and Option-Lok®. Latex-Free.	104" (264 cm)	15	17.0	0.090	48
						
Primary I.V. Sets – Micro bore						
16000-01	LifeShield® Symbiq™ Pump Set, Microdrip® – with non-DEHP tubing, convertible pin, piggyback with backcheck valve, 2 CLAVE® Y-sites, CAIR® clamp, slide clamp, and Option-Lok®. Latex-Free.	106" (269 cm)	60	9.6	0.054	48
						

List No.	Product Description	Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size	List No.	Product Description	Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
Primary I.V. Sets – Microbore (cont.)							Burette Sets – Microdrip						
16002-01	LifeShield® Symbiq™ Pump Set – with yellow striped tubing, convertible pin, slide clamp, distal microbore tubing, and Option-Lok®. Latex-Free.	108.5" (276 cm)	15	8.8	0.054	48	16004-01	LifeShield® Symbiq™ Pump 50 mL Burette Set, Microdrip®, SoluSet® – with non-DEHP tubing, convertible pin, 3 slide clamps, bacterial retentive air filter with air filter slide clamp, CLAVE® port, burette, CAIR® clamp, 2 CLAVE® Y-sites, distal microbore tubing, and Option-Lok®. Latex-Free.	104" (264 cm)	60	8.1	0.054	24
													
16024-01	LifeShield® Symbiq™ Pump Set – with orange polyethylene-lined light resistant tubing, convertible pin, slide clamp, CLAVE® Y-site, distal microbore tubing, and Option-Lok®. Latex-Free.	105" (267 cm)	15	8.3	0.043	48	16009-01	LifeShield® Symbiq™ Pump 50 mL Burette Set, Microdrip®, SoluSet® – with non-DEHP tubing, convertible piercing pin, 2 slide clamps, bacterial retentive air filter with clamp, CLAVE® port, burette, CAIR® clamp, 0.2 micron filter, CLAVE® Y-site, distal microbore tubing, and Option-Lok®. Latex-Free.	100" (254 cm)	60	10.7	0.054	24
													
Burette Sets – Macro drip													
16003-01	LifeShield® Symbiq™ Pump 150 mL Burette Set, SoluSet® – with non-DEHP tubing, convertible pin, 3 slide clamps, bacterial retentive air filter with slide clamp, CLAVE® port, burette, CAIR® clamp, 2 CLAVE® Y-sites, 0.2 micron filter, and Option-Lok®. Latex-Free.	117" (297 cm)	15	17.0	0.090	24	16014-01	LifeShield® Symbiq™ Pump 150 mL Burette Set, Microdrip®, SoluSet® – with non-DEHP tubing, convertible pin, 3 slide clamps, bacterial retentive air filter with slide clamp, upper CLAVE® port, burette, CAIR® clamp, 2 CLAVE® Y-sites, distal microbore tubing, and Option-Lok®. Latex-Free.	124.5" (316 cm)	60	9.5	0.054	24
													
16017-01	LifeShield® Symbiq™ Pump 150 mL Burette Set, SoluSet® – with non-DEHP tubing, convertible pin, slide clamp, bacterial retentive air filter with slide clamp, CLAVE® port, burette, CAIR® clamp, 2 CLAVE® Y-sites, and Option-Lok®. Latex-Free.	115" (292 cm)	15	15.0	0.090	24	16018-01	LifeShield® Symbiq™ Pump 150 mL Burette Set, Microdrip®, SoluSet® – with non-DEHP tubing, convertible pin, 4 slide clamps, bacterial retentive air filter with slide clamp, CLAVE® port, burette, CAIR® clamp, 3 CLAVE® Y-sites, 0.2 micron filter, distal microbore tubing and Option-Lok®. Latex-Free.	138" (351 cm)	60	10.0	0.054	24
													

List No.	Product Description	Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
Blood Administration Sets						
16006-01	LifeShield® Symbiq™ Pump Blood Set – with nonvented pin, 210 micron blood filter, CAIR® clamp, CLAVE® Y-site, and Option-Lok®. Latex-Free.	110" (279 cm)	10	66.0	0.090	24
						
16007-01	LifeShield® Hemoset® Symbiq™ Pump 100 mL Burette Set – with non-DEHP tubing, piercing pin, slide clamp, bacterial retentive air filter with slide clamp, CLAVE® port, burette with 170 Micron Filter, CAIR® clamp, and Option-Lok®. Latex-Free.	104" (264 cm)	15	14.0	0.090	24
						
16020-01	LifeShield® Symbiq™ Pump Hema® Y-Type Blood Set – with 2 nonvented pins, 3 CAIR® clamps, 210 micron blood filter, CLAVE® Y-site, and Option-Lok®. Latex-Free.	122" (310 cm)	10	66.0	0.090	24
						
Enteral Sets						
16015-01	Symbiq™ Pump Enteral Set – with 40 mm screw cap with non-DEHP tubing, bacterial retentive air filter, CAIR® clamp, and non-I.V. enteral catheter adapter. Latex-Free.	99" (251 cm)	15	12.0	0.090	24
						
Nitroglycerin Sets						
16016-01	Symbiq™ Enteral Pump Set – with non-DEHP tubing, filling port, integral feeding container, CAIR® clamp, and non-I.V. enteral catheter adapter. Latex-Free.	102" (259 cm)	15	11.0	0.090	24
						
16010-01	Symbiq™ Pump Nitroglycerin Set – with polyethylene-lined tubing, convertible pin, slide clamp, and Option-Lok®. Latex-Free.	108.5" (276 cm)	15	16.4	0.100	48
						
16025-01	Symbiq™ Pump Nitroglycerin Set, Microdrip® – with polyethylene-lined tubing, convertible pin, slide clamp, and Option-Lok®. Latex-Free.	108" (274 cm)	60	17.2	0.100	48
						
Syringe Sets						
16023-01	LifeShield® Symbiq™ Pump Syringe Microbore Set – with non-DEHP tubing, upper CLAVE® port, syringe holder, slide clamp, CLAVE® Y-site, and Option-Lok®. Latex-Free.	83" (211 cm)	-	3.5	0.054	56
						

Symbiq™ Infusion System with Hospira MedNet® Software



NON-DEHP				NON-DEHP in FLUID PATH		DEHP in FLUID PATH	
16000-01	16003-01	16004-01	16005-01	16010-01	16013-01	16001-01	16002-01
16007-01	16009-01	16011-01	16014-01	16024-01	16025-01	16006-01	16008-01
16015-01	16017-01	16018-01	16019-01			16020-01	
16021-01	16022-01	16023-01					



Advancing Wellness™

www.hospira.com

©Hospira, Inc. 275 North Field Drive, Lake Forest, IL 60045 P06-0589-Nov, 06