

IHE INTEGRATION STATEMENT

IHE INTEGRATION STATEMENT¹			
VENDOR	PRODUCT NAME	VERSION	DATE
HOSPIRA	HOSPIRA MEDNET™	6.1	JULY 2015
THIS PRODUCT IMPLEMENTS ALL TRANSACTIONS REQUIRED IN THE IHE TECHNICAL FRAMEWORK TO SUPPORT THE IHE INTEGRATION PROFILES, ACTORS AND OPTIONS LISTED BELOW:			
INTEGRATED PROFILES IMPLEMENTED	ACTORS IMPLEMENTED	OPTIONS IMPLEMENTED	
CT – CONSISTENT TIME	TIME CLIENT		
DEC – DEVICE ENTERPRISE COMMUNICATION	DEVICE OBSERVATION REPORTER		
PIV – POINT-OF-CARE INFUSION VERIFICATION	INFUSION ORDER CONSUMER		
ACM – ALERT COMMUNICATION MANAGEMENT	ALERT REPORTER		
IPEC – INFUSION PUMP EVENT COMMUNICATION	DEVICE OBSERVATION REPORTER		
INTERNET ADDRESS FOR VENDOR'S IHE INFORMATION:		http://www.ivclinicalintegration.com/news-resources	
LINKS TO STANDARDS CONFORMANCE STATEMENTS FOR IMPLEMENTATION			
HL7	NOT APPLICABLE		
DICOM	NOT APPLICABLE		
LINKS TO GENERAL INFORMATION ON IHE			
IN NORTH AMERICA:	IN EUROPE:	IN JAPAN:	
www.ihe.net	www.ihe-europe.net	www.jira-net.or.jp/ihe-j	

1. Hospira, Inc. Data on File. 2015 IHE Connectathon Result Matrix. March 2015.

